



For Immediate Release

Metrohm 885 Compact Oven Sample Changer Wins 2011 Laboratory Equipment Reader's Choice Award

Riverview, FL, February 2011...

The readers of Laboratory Equipment Magazine have chosen the 885 Compact Oven Sample Changer as the winner in the Moisture and Particle Analyzer Category. Our newest KF Oven makes Karl Fischer titrations easier than ever before. It's perfect for any sample that you want to analyze for water content—from low ppm to 100%! The model 885 has the following exclusive features:

Metrohm's unique headspace KF technology

Often imitated but never duplicated, our innovative technology for analyzing moisture content via the KF oven technique uses headspace vials as the oven chamber! No more messy tube ovens or drying balances!

Small footprint – plenty of room for samples

It's right there in the name: 885 Compact Oven Sample Changer. Don't judge it by its size, the 885 has the footprint of a dictionary, but you can fill the rack with 18 samples, standards and QC check samples—it handles them all!

Great flexibility

No matter if you perform coulometric or volumetric Karl Fischer titrations, the 885 Compact Oven Sample Changer can be operated with all stand-alone KF titrators.

Easy to use

A comfortable and user-friendly keypad controls the 885 Compact Oven Sample Changer: just input the number of samples, your oven temperature and the gas flow – that's all!

Publishing Contact:

Christy Hurlles
866-METROHM (toll-free)
813-316-4900 (fax)
churlles@metrohmusa.com

Press relations for:

Metrohm USA, Inc.
Metrohm Canada, Inc.
Applikon USA, a Metrohm company

9702.N1.1009

Metrohm USA, Inc. • Tel 800-727-6768 • Fax 813-316-4900 • Email info@metrohmusa.com • www.metrohmusa.com