

# Make eluent “to go”



With premixed concentrates that are convenient—and affordable.

**Because \$2.40 and 2 seconds are all it takes to make 1L of certifiable, properly mixed eluent,** join the many labs already making the switch from “home brewing.” Use the time you would otherwise spend weighing, pipetting and mixing to take on more important work in your lab. And rest assured that your runs—as well as your results—are always optimized.



SNIPS eluent concentrates are Certified & NIST-traceable. Turn over for more details.

 **Metrohm**  
USA Inc.

# Metrohm® Eluent Concentrate SNIPS

Each stay-fresh Teflon®-sealed tube contains 10 mL of premeasured, certified eluent concentrate. We call them “SNIPS” because you simply “snip off” an end and pour into enough DI water to make 1L of perfectly mixed eluent. For added convenience, SNIPS are available in packs of 25 and 50.



## SNIPS provide the following benefits:

- Certified pre-measured concentrates ensure accurate runs and consistent results
- Sealed, inert Teflon packaging provides long (18-month) shelf life
- Eliminate time spent weighing and mixing to get to more important work
- No more run errors and poor results due to incorrect weighing and pipetting
- Conveniently available for A Supp 5 and A Supp 7 eluents in packs of 25 or 50

## Ordering Information:

Catalog No.	Description	List USD
M.920.0562	A Supp 5 Eluent, 100x concentrate, 10mL SNIPS, pack of 25	\$60.00
M.920.0562k	A Supp 5 Eluent, 100x concentrate, 10mL SNIPS, pack of 50	\$120.00
M.920.0567	A Supp 7 Eluent, 100x concentrate, 10mL SNIPS, pack of 25	\$60.00
M.920.0567k	A Supp 7 Eluent, 100x concentrate, 10mL SNIPS, pack of 50	\$120.00

866-METROHM (638-7646)  
Web: [www.metrohmusa.com](http://www.metrohmusa.com)  
Email: [info@metrohmusa.com](mailto:info@metrohmusa.com)

 **Metrohm**  
USA Inc.